



Micro+™ Smokerlyzer®

Helping people to stop smoking, one breath at a time.




One CO monitor for all your testing needs



FREE
COdata+™
software

Ideal for in-depth stop smoking advice delivered by:

- Stop smoking services
- General practice
- Schools
- Pharmacies
- Health visitors
- Midwives
- Researchers



Automatic %FCOHb conversion



Create and save patient details



For ambient levels of CO

Smokerlyzer®
www.bedfont.com

Features and benefits

- Brand new easy to use interface
- Create and store patient details
- SteriTouch® technology for improve infection control
- Large touch screen
- Automatic conversion to %FCOHb
- Familiar green, amber and red traffic light system
- Ambient levels of CO

Technical specification

Concentration range	0-500ppm
Display	Full colour touchscreen
Detection principle	Electrochemical sensor
Repeatability	<±5%
Accuracy	2ppm/±5%*
Power	3 x AA (LR6 or equivalent) - up to 1000 minutes 1 x CR2032 Lithium coin cell
T₉₀ response time	<30 seconds
Operating temperature	15-40°C
Storage/transport temperature	0-50°C
Operating/storage/transport pressure	Atmospheric ±10%
Operating humidity	15-90% non-condensing
Storage/transport humidity	0-95%
Sensor operating life	5 years
Sensor sensitivity	1ppm
Sensor drift	<5% per annum
Dimensions	Approx. 37 x 77 x 140mm
Weight	Approx. 215g (including batteries)
Materials	Case: polycarbonate/ABS blend with SteriTouch® anti-microbial additive D-piece™: polypropylene SteriBreath™: polypropylene
H₂ cross interference	≤6%

*Whichever is greater



British technology

Bedfont® Scientific Ltd

Station Road, Harrietsham, Maidstone
Kent, ME17 1JA, England
Tel: +44 (0)1622 851122, Fax: +44 (0)1622 854860
Email: ask@bedfont.com
www.bedfont.com

© Bedfont® Scientific Limited 2017

Issue 8 - April 2017 - Part No: MKT375
Bedfont® Scientific Limited reserve the right to change or update this literature without prior notice.
Registered in: England and Wales. Registered No: 1289798



DUOTEC AS
HÅLEINSTRUMENTER TIL MILJØKONTROL

DUOTEC A/S
Herstedøstervej 19 · 2600 Glostrup
Tlf. 4345 9188 · Fax 4343 2330
e-mail: info@duotec.dk
www.duotec.dk

bedfont®
est. 1976